

MIXPAC™

T-Mixer™ with Intra Oral Tip for 18 to 75 mL Systems



Mixing tips Indications

- Impression Taking
- Bite Registrations
- Temporary Crown & Bridge
- Permanent Cementation
- Post Cementation & Core Build-ups
- Relining

Now available in additional Candy Colors™

IOT 406-70

Ratio	MOQ Box	Outer Ø (mm)	Inner Ø (mm)
1:1 / 2:1	70'000	1.4	0.6

MBT 4.2-12-D

Ratio	MOQ Box	Inner Ø (mm)	Elements	L (mm)	Vol. (mL)
1:1 / 2:1	8'000	4.2	12	50.1	0.59



IOT 506-70

Ratio	MOQ Box	Outer Ø (mm)	Inner Ø (mm)
1:1 / 2:1	70'000	1.4	0.6

MBT 5.4-12-D

Ratio	MOQ Box	Inner Ø (mm)	Elements	L (mm)	Vol. (mL)
1:1 / 2:1	6'000	5.4	12	59.4	1.02




IOT 506-70

Ratio	MOQ Box	Outer Ø (mm)	Inner Ø (mm)
1:1 / 2:1	70'000	1.4	0.6

MBT 6.5-12-D

Ratio	MOQ Box	Inner Ø (mm)	Elements	L (mm)	Vol. (mL)
1:1 / 2:1	5'000	6.5	12	67.6	1.60



MIXPAC owns valuable trademark protection for the shape, unique looking X marking and color of its products.



All MIXPAC Candy Colors™ are available upon request for Intra Oral Tips; and mixing tips in any diameter, length and configuration for diverse material types and applications. The Candy Colors stand for our Swiss Quality.

medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland
P +41 81 414 70 00 | www.medmix.swiss



DENTAL

Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, T-Mixer, and Candy Colors are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023